

Catalogue Number	Product	Order number / Unit
1130	ACY-FDG Reference standard for byproduct of [¹⁸F]FDG synthesis (2-[¹⁸F]Fluoro-2-deoxy-D-glucose) Molar Mass: 350.29 $C_{14}H_{19}FO_9$ [141395-48-4] Yellowish solid packaged in dark glass screw cap vials. Purity: > 95 % Certificates: CoA; ¹ H and ¹⁹ F NMR spectra Chemical Name: CA index name: D-Glucopyranose, 2-deoxy-2-fluoro-, tetraacetate Synonymes: 2-Fluoro-2-deoxy-glucose tetraacetate Literature: 1. Kovac P. et al. Short synthesis of 2-Deoxy-2-fluoro-D-glucose. Carbohydrate Research 1986, 153, 168-170.	1130.0010: 10 mg per vial 1130.0020: 20 mg per vial Please inquire for customized filling and bulk quantities. 