

Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

1131**1,2,3,4-tetra-O-acetyl-beta-D-glucopyranose****TLC Standard for Testing of Radiochemical Purity of [¹⁸F]FDG according to Ph. Eur. 6.2, 1325****Molar Mass:** 348.3C₁₄H₂₀O₁₀

[13100-46-4]

Colourless solid packaged in dark glass crimp cap vials.

Purity: > 95 %**Certificates:**CoA; ¹H NMR and ¹³C NMR spectra**Chemical Name:**

CA index name: beta-D-Glucopyranose, 1,2,3,4-tetraacetate

Synonymes:

Glucopyranose, 1,2,3,4-tetraacetate, β-D-; 1,2,3,4-Tetra-O-acetyl-β-D-glucopyranose; 1,2,3,4-Tetra-O-acetyl-β-D-glucose; 6-OH-TAG

Literature:

European Pharmacopoeia 6.2, 3678-3680.

1131.0030: 30 mg per vial
Please inquire for customized filling and bulk quantities.