


Catalogue Number	Product	Order number / Unit
1140	Talose Triflate Precursor for [¹⁸F]FDGal (2-[¹⁸F]Fluoro-2-deoxy-D-galactose) Molar Mass: 480.37 $C_{15}H_{19}F_3O_{12}S$ CAS-RN not yet assigned Colourless or nearly colourless crystals packaged in dark glass crimp cap vials. Purity: > 95 % Certificates: CoA; ¹ H and ¹⁹ F NMR spectra Chemical Name: 1,3,4,6-Tetra-O-acetyl-2-O-trifluoromethanesulfonyl-β-D-talopyranose Synonymes: β-D-talopyranose, 1,3,4,6-tetraacetate 2-(trifluoromethanesulfonate) Literature: 1. Fukuda H. et al. (F-18)-2-Deoxy-2-Fluoro-D-Galactose: A Potential Tracer for the Assessment of Sugar Metabolism in the Liver by Positron Emission Tomography. Cyric Annual Report 1982, 188-195.	1140.0010: 10 mg per vial 1140.0025: 25 mg per vial Please inquire for customized filling and bulk quantities. 

date of product catalogue issue: 22 December 2010