
Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

1231 **5'-O-Benzoyl-2,3'-anhydrothymidine (GMP)**

1231.0010: 10 mg per vial
Please inquire for customized filling and bulk quantities.

Precursor for [¹⁸F]FLT (3'-Deoxy-3'-[¹⁸F]fluorothymidine)

Manufactured according to GMP requirements for APIs for use in clinical trials (ICH Q7 chapter 19)

Molar Mass: 328.32

C₁₇H₁₆N₂O₅

[70838-44-7]

White solid packaged in dark glass crimp cap vials.

Purity: > 95 %

Certificates:

CoA with ¹H, ¹³C NMR and IR spectra (identity); HPLC (purity); HS-GC (residual solvents)

Chemical Name:

For chemical information please refer to product number 1230.

Synonymes:

Please refer to product number 1230.

Literature:

For literature references please refer to product number 1230.

