

Catalogue Number	Product	Order number / Unit
1241	<p><b>3-N-Boc-5'-O-dimethoxytrityl-3'-O-nosyl-thymidine (GMP)</b></p> <p><b>Precursor for [<sup>18</sup>F]FLT (3'-Deoxy-3'-[<sup>18</sup>F]fluorothymidine)</b></p> <p>Manufactured according to GMP requirements for APIs for use in clinical trials (ICH Q7 chapter 19)</p> <p><b>Molar Mass:</b> 829.87</p> <p>C<sub>42</sub>H<sub>43</sub>N<sub>3</sub>O<sub>13</sub>S</p> <p>[444717-23-1]</p> <p>Yellowish amorphous solid packaged in clear glass vials. (10 ml headspace)</p> <p><b>Purity:</b> ≥ 97 %</p> <p><b>Certificates:</b> CoA with <sup>1</sup>H, <sup>13</sup>C NMR and IR spectra (identity); HPLC (purity); HS-GC (residual solvents); Content Silver (AAS); tested for Bioburden and Bacterial endotoxins</p> <p><b>Chemical Name:</b> For chemical information please refer to product number 1240.</p> <p><b>Synonymes:</b> Please refer to product number 1241.</p> <p><b>Literature:</b> For literature references please refer to product number 1240.</p>	<p>1241.0025: 25 mg per vial 1241.0040: 40 mg per vial Please inquire for customized filling and bulk quantities.</p> 