


Catalogue Number	Product	Order number / Unit
1380	6-Fluorodopamine Reference standard for 6-[¹⁸F]Fluorodopamine Molar Mass: 171.17 $C_8H_{10}FNO_2$ [71144-39-3] Colourless solid packaged in dark glass screw cap vials. Purity: > 95 % Certificates: CoA; ¹ H and ¹⁹ F NMR spectra Chemical Name: CA index name: 1,2-Benzenediol, 4-(2-aminoethyl)-5-fluoro- Synonymes: 6-Fluorodopamine; 4-(2-Aminoethyl)-5-fluorobenzene-1,2-diol; 4-Fluoro-5-(2-aminoethyl)benzene-1,2-diol Literature: 1. Ding et al. Synthesis of High Specific Activity 6-[¹⁸ F]Fluorodopamine for Positron Emission Tomography Studies of Sympathetic Nervous Tissue. J. Med. Chem. 1991, 34. 861-863.	1380.0005: 5 mg per vial 1380.0010: 10 mg per vial Please inquire for customized filling and bulk quantities. 

date of product catalogue issue: 05 April 2012