


Catalogue Number	Product	Order number / Unit
1591	Fluoroethyl-spiperone Reference standard for 3-(2'-[¹⁸F]Fluorethyl) spiperone Molar Mass: 441.51 $C_{25}H_{29}F_2N_3O_2$ [106114-42-5] Colourless solid packaged in dark glass screw cap vials. Purity: > 95 % Certificates: CoA; ¹ H and ¹⁹ F NMR spectra Chemical Name: CA index name: 3-(2-fluoroethyl)-8-[4-(4-fluorophenyl)-4-oxobutyl]-1-phenyl-1,3,8-triazaspiro[4.5]decan-4-one Synonymes: FESP, Fluoroethylspiperone Literature: <ol style="list-style-type: none">1. Moerlein S.M. et al. Effect of Lipophilicity on the in Vivo Localization of Radiolabelled spiperone Analogues Int. J. Nucl. Med. Biol. 1985, 12, 353-356.2. Chi D.Y. et al. Synthesis of no carrier added N-(3-[¹⁸F]fluoroalkyl)-spiperone derivatives. Appl. Radiat. Isot. 1986, 37, 1173-1180.3. Coenen H.H. et al. 3-N-(2-[¹⁸F]Fluorethyl)-spiperone: A novel ligand for cerebral dopamine receptor studies with PET Life Sciences 1987, 40, 81-88.4. Satyamurthy N. et al. 3-(2-[¹⁸F]Fluorethyl)spiperone, a Potent Dopamine Antagonist: Synthesis, Structural Analysis and In-vivo Utilization in Humans Int. J. Rad. Appl. Instrum. A 1990, 41, 113-129.5. Hamacher K. et al. Remote controlled one-step production of ¹⁸F-labeled butyrophenone neuroleptics exemplified by the synthesis of n.c.a. [¹⁸F]N-methylspiperone. Appl. Radiat. Isot. 1995, 46, 911-916.	1591.0010: 10 mg per vial Please inquire for customized filling and bulk quantities. 

date of product catalogue issue: 05 April 2012