

Catalogue Number	Product	Order number / Unit
1655	<b>FEPPA Standard</b> Reference standard for [ <sup>18</sup> F]-FEPPA ([ <sup>18</sup> F]-N-(2-(2-Fluoroethoxy)benzyl)-N-(4-phenoxy-3-pyridin-3-yl)acetamide) <b>Molar Mass:</b> 380.41 C <sub>22</sub> H <sub>21</sub> FN <sub>2</sub> O <sub>3</sub> [1095267-90-5] Colourless to yellowish oil packaged in dark glass screw cap vials. <b>Purity:</b> > 95 % <b>Certificates:</b> CoA; <sup>1</sup> H and <sup>19</sup> F NMR spectra <b>Chemical Name:</b> CA index name: Acetamide, N-[[2-(2-fluoroethoxy)phenyl]methyl]-N-(4-phenoxy-3-pyridinyl)- <b>Synonymes:</b> N-Acetyl-N-(2-fluoroethoxybenzyl)-2-phenoxy-5-pyridinamine <b>Literature:</b> 1. Wilson A.A. et al. Radiosynthesis and initial evaluation of [ <sup>18</sup> F]-FEPPA for PET imaging of peripheral benzodiazepine receptors. Nucl. Med. Biol. 2008, 35, 305-314.	1655.0010: 10 mg per vial Please inquire for customized filling and bulk quantities. 