

Catalogue Number	Product	Order number / Unit
1721	Tosyloxyethylflumazenil Precursor for [¹⁸F]FFMZ (2'-[¹⁸F]Fluorethylflumazenil) Molar Mass: 473.47 C₂₂H₂₀FN₃O₆S [676437-17-5] Colourless solid packaged in dark glass crimp cap vials. Purity: > 95 % Certificates: CoA; ¹ H and ¹⁹ F NMR spectra Chemical Name: CA index name: 4H-Imidazo[1,5-a][1,4]benzodiazepine-3-carboxylic acid, 8-fluoro-5,6-dihydro-5-methyl-6-oxo-, 2-[[[(4-methylphenyl)sulfonyl]oxy]ethyl] ester Synonymes: TFMZ; 8-Fluoro-5-methyl-6-oxo-5,6-dihydro-4H-2,5,10b-triaza-benzo[e]azulene-3-carboxylic acid 2-(toluene-4-sulfonyloxy)-ethyl ester; 2'-Tosyloxyflumazenil Literature: 1. Yoon Y.H. et al. Novel one-pot one-step synthesis of 2'-[¹⁸ F]fluoroflumazenil (FFMZ) for benzodiazepine receptor imaging. Nucl. Med. Biol. 2003, 30, 521-527.	1721.0005: 5 mg per vial 1721.0010: 10 mg per vial Please inquire for customized filling and bulk quantities. 