


Catalogue Number	Product	Order number / Unit
1910	16alpha-Fluoroestradiol Reference standard for 16alpha-[¹⁸F] Fluoroestradiol Molar Mass: 290.37 $C_{18}H_{23}FO_2$ [92817-10-2] Colourless solid packaged in dark glass screw cap vials. Purity: > 90 % Certificates: CoA; ¹ H and ¹⁹ F NMR spectra Chemical Name: CA index name: Estra-1,3,5(10)-triene-3,17-diol, 16-fluoro-, (16alpha,17beta) Synonymes: 16alpha-Fluoro-13beta-methyl-1,3,5(10)-gonatriene-3,17beta-diol; 16alpha-Fluoro-17beta-estradiol; FES	1910.0001: 1 mg per vial 1910.0002: 2 mg per vial 1910.0010: 10 mg per vial Please inquire for customized filling and bulk quantities.



date of product catalogue issue: 05 April 2012