

Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

2510 (S)-Desisopropylcarazolol

Precursor for (S)-[¹¹C]Carazolol

Molar Mass: 256.3

C₁₅H₁₆N₂O₂

[143412-40-2]

Colourless or nearly colourless solid packaged in dark glass crimp cap vials.

Purity: > 95 %

Certificates:

CoA; ¹H and ¹³C NMR spectra

Chemical Name:

CA index name: 2-Propanol, 1-amino-3-(9H-carbazol-4-yloxy)- (S)


Synonymes:

(S)-1-amino-3-(9H-carbazol-4-yloxy)-2-propanol

Literature:

1. Berridge M.S. et al. Preparation and in Vivo Binding of [¹¹C]Carazolol, a Radiotracer for the Beta-adrenergic Receptor. Nucl. Med. Biol. 1992, 19, 563-569.
2. Zheng L. et al. Synthesis, Binding Properties, and ¹⁸F-Labeling of Fluorocarazolol, a High-Affinity beta-Adrenergic Receptor Antagonist. J. Med. Chem. 1994, 37, 3219-3230.
3. Elsinga P.H. et al. (S,S)- and (S,R)-1'-[¹⁸F]Fluorocarazolol, Ligands for the Visualization of Pulmonary beta-Adrenergic Receptors with PET. Nucl. Med. Biol. 1996, 23, 159-167.

2510.0002: 2 mg per vial
Please inquire for customized filling and bulk quantities.



date of product catalogue issue: 05 April 2012