
Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

2952

Nitro-Xeloda

Precursor for [¹⁸F]Xeloda

Molar Mass: 470.43

C₁₉H₂₆N₄O₁₀

[865474-03-9]

Colourless to yellowish solid packaged in dark glass screw cap vials.

Purity: > 95 %

Certificates:

CoA; ¹H and ¹³C NMR spectra

Chemical Name:

CA index name: Cytidine, 5'-deoxy-5-nitro-N-[(pentoxy)carbonyl]-, 2',3'-diacetate


Synonyms:

2',3'-Di-O-acetyl-5'-deoxy-5-nitro-N4-(pentyloxycarbonyl)cytidine; Nitroxeloda

Literature:

1. Xiangshu F. et al. Synthesis of [¹⁸F]Xeloda as a novel potential PET radiotracer for imaging enzymes in cancers. Nucl. Med. Biol. 2004, 31, 1033-1041.

2952.0010: 10 mg per vial
Please inquire for customized filling and bulk quantities.



date of product catalogue issue: 05 April 2012