

Catalogue Number	Product	Order number / Unit
3180	2-Iodo-alpha-methyl-p-tyrosine Reference standard for 2-[¹²³I]Iodo-α-methyl-p-L-tyrosine Molar Mass: 321.11 $C_{10}H_{12}INO_3$ CAS-RN not yet assigned Off-white solid packaged in dark glass screw cap vials. Purity: 95 % (HPLC) Certificates: CoA; ¹ H NMR spectrum Chemical Name: L-Tyrosine, 2-iodo-α-methyl- Synonymes: Tyrosine, 2-iodo-α-methyl-, L- ; L-2-Iodo-α-methyl-tyrosine ; o-Iodo-α-methyl-L-tyrosine, 2-IAMT Literature: no literature reference available	3180.0010: 10 mg per vial Please inquire for customized filling and bulk quantities. 