

Catalogue Number	Product	Order number / Unit
<b>4211</b>	<b>Altropane hydrochloride</b> <b>Reference standard for [<sup>11</sup>C]Altropane</b> <b>Reference standard for [<sup>125</sup>I]Altropane</b>  <b>Molar Mass:</b> 465.73 C <sub>18</sub> H <sub>21</sub> FINO <sub>2</sub> · HCl  CAS-RN not yet assigned Altropane free amine: [180468-34-2] Colourless solid packaged in dark glass screw cap vials.  <b>Purity:</b> > 95 % <b>Certificates:</b> CoA; <sup>1</sup> H and <sup>19</sup> F NMR spectra <b>Chemical Name:</b> 8-Azabicyclo[3.2.1]octane-2-carboxylic acid, 3-(4-fluorophenyl)-8-[(2E)-3-iodo-2-propenyl]-, methyl ester, (1R,2S,3S,5S)-, hydrochloride salt <b>Synonymes:</b> (E)-N-(1-Iodoprop-1-en-3-yl)-3-β-(4-fluorophenyl)-nortropane-2-β-carboxylic acid methyl ester, hydrochloride; (E)-N-(3-Iodoprop-2-enyl)-2β-carbomethoxy-3β-(4-fluorophenyl)-nortropane hydrochloride; IACFT HCl <b>Literature:</b> For literature references please refer to Prod. Nr. 4210	4211.0001: 1 mg per vial 4211.0010: 10 mg per vial Please inquire for customized filling and bulk quantities.  