

Catalogue  
Number

Product

Order number / Unit

**4231**

**Benztropine Mesylate**

Reference standard for [<sup>11</sup>C]Benztropine

**Molar Mass:** 403.53

C<sub>21</sub>H<sub>25</sub>NO · CH<sub>4</sub>O<sub>3</sub>S

[132-17-2]

Colourless solid packaged in dark glass screw cap vials.

**Purity:** > 95 %

**Certificates:**

CoA; <sup>1</sup>H spectrum

**Chemical Name:**

3-(benzhydryloxy)-8-methyl-8-azabicyclo[3.2.1]octane methanesulfonate


**Synonymes:**

8-azabicyclo[3.2.1]octane, 3-(diphenylmethoxy)-, endo, methanesulfonate;

**Literature:**

1. Agoston G.E. et al. Novel N-Substituted 3 $\alpha$ -[Bis(4'-fluorophenyl)methoxy]tropane Analogues: Selective Ligands for the Dopamine Transporter. J. Med. Chem. 1997, 40, 4329-4339

Please inquire for customized filling and bulk quantities.



date of product catalogue issue: 05 April 2012