

---

Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

---

**4301****MASB (GMP)****Precursor for [<sup>11</sup>C]DASB  
([<sup>11</sup>C]N,N-Dimethyl-2-(2-amino-4-cyanophenylthio)  
benzylamine)**Manufactured according to GMP requirements for APIs  
for use in clinical trials (ICH Q7 chapter 19)**Molar Mass:** 269.37C<sub>15</sub>H<sub>15</sub>N<sub>3</sub>S

[296774-10-2]

Off-white to brownish solid packaged in dark glass  
screw cap vials**Purity:** > 98 %**Certificates:**CoA with <sup>1</sup>H, <sup>13</sup>C NMR, IR spectra, and melting point  
(identity); HPLC (purity); HS-GC (residual solvents)**Chemical Name:**

For chemical information please refer to product number 4300.

**Synonymes:**

Please refer to product number 4300.

**Literature:**

For literature references please refer to product number 4300.

4301.0050: 50 mg per vial  
Please inquire for customized  
filling and bulk quantities.