
Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

4311 DASB (GMP)

Reference standard for [¹¹C]DASB
([¹¹C]N,N-Dimethyl-2-(2-amino-4-cyanophenylthio)
benzylamine)

Manufactured according to GMP requirements for APIs
for use in clinical trials (ICH Q7 chapter 19)

Molar Mass: 283.39

C₁₆H₁₇N₃S

[627490-01-1]

Yellowish solid packaged in dark glass screw cap vials

Purity: > 98 %

Certificates:

CoA with ¹H, ¹³C NMR, IR spectra, and melting point
(identity); HPLC (purity); HS-GC (residual solvents)

Chemical Name:

For chemical information please refer to product number 4310.

Synonymes:

Please refer to product number 4310.

Literature:

For literature references please refer to product number 4310.

4311.0050: 50 mg per vial
Please inquire for customized
filling and bulk quantities.

