

---

Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

---

**4397****FPyKYNE****Reference standard for [<sup>18</sup>F]FPyKYNE  
2-[<sup>18</sup>F]Fluoro-3-pent-4-yn-1-yloxy pyridine****Molar Mass:** 179.19C<sub>10</sub>H<sub>10</sub>FNO

[1161009-63-7]

Yellow solid packaged in dark glass screw cap vials.

**Purity:** > 95 %**Certificates:**CoA; <sup>1</sup>H NMR spectrum**Chemical Name:**

CA index name: Pyridine, 2-fluoro-3-(4-pentyn-1-yloxy)-

**Synonymes:**

FPyKYNE; 2-Fluoro-3-pent-4-ynyloxy-pyridine

**Literature:**

1. Dollé F. et al. [<sup>18</sup>F]FPyKYNE, a fluoropyridine-based alkyne reagent designed for the fluorine-18 labelling of macromolecules using click chemistry. J. Labelled Compd. Radiopharm. 2008, 51, 336-342.

**4397.0010: 10 mg per vial  
Please inquire for customized  
filling and bulk quantities.**