

Catalogue Number	Product	Order number / Unit
<b>4406</b>	<b>(2S,3S)-Methylreboxetine</b> <b>Reference standard for [<sup>11</sup>C]MeNER</b> <b>( [<sup>11</sup>C]Methyl norethyl reboxetine)</b>  <b>Molar Mass:</b> 299.36 $C_{18}H_{21}NO_3$ [105018-14-2]  Colourless to white solid packaged in dark glass crimp cap vials.  <b>Purity:</b> > 95 %  <b>Certificates:</b> CoA; <sup>1</sup> H NMR spectrum  <b>Chemical Name:</b> CA index names: Morpholine, 2-[(S)-(2-methoxyphenoxy)phenylmethyl]-, (2S)-; Morpholine, 2-[(2-methoxyphenoxy)phenylmethyl]-, [S-(R*,R*)]-  <b>Synonymes:</b> (S,S)-2-( $\alpha$ -(2-Methoxyphenoxy)benzyl)morpholine; Methyl norethyl reboxetine; MeNER  <b>Literature:</b> 1. Schou M. et al. Specific in vivo binding to the norepinephrine transporter demonstrated with the PET radioligand, (S,S)-[ <sup>11</sup> C]MeNER. Nucl. Med. Biol. 2003, 30, 707-714.	4406.0010: 10 mg per vial 4406.0100: 100 mg per vial Please inquire for customized filling and bulk quantities.

