
Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

5161 3-(2-Benzoxazolyl)-7-hydroxycoumarin

Precursor for 3-(2-Benzoxazolyl)-7-[¹¹C] methoxycoumarin

Molar Mass: 279.25

C₁₆H₉NO₄

[64887-40-7]

Yellow to orange solid packaged in dark glass screw cap vials.

Purity: > 95 %

Certificates:

CoA; ¹H NMR spectrum

Chemical Name:

CA index name: 2H-1-Benzopyran-2-one, 3-(2-benzoxazolyl)-7-hydroxy-

Synonymes:

3-(2-Benzoxazolyl)umbelliferone; 3-(2-Benzoxazolyl)-7-hydroxycoumarin; 3-Benzoxazol-2-yl-7-hydroxy-chromen-2-one

Literature:

no literature reference available

5161.0010: 10 mg per vial
5161.0100: 100 mg per vial
Please inquire for customized filling and bulk quantities.

