


Catalogue Number	Product	Order number / Unit
5610	(+)-Norverapamil Precursor for (+)-[¹¹C]Verapamil Molar Mass: 440.58 $C_{26}H_{36}N_2O_4$ [123932-43-4] Yellowish to brownish oil packaged in dark glass crimp cap vials. Purity: > 95 % Certificates: CoA; ¹ H and ¹³ C NMR spectra Chemical Name: CA index name: Benzeneacetonitrile, α -[3-[[2-(3,4-dimethoxyphenyl)-ethyl]amino]propyl]-3,4-dimethoxy- α -(1-methylethyl)-, (α R)- Synonyms: (α R)- α -[3-[[2-(3,4-Dimethoxyphenyl)ethyl]amino]propyl]-3,4-dimethoxy- α -(1-methylethyl)-benzeneacetonitrile; (+)-Norverapamil; R-Norverapamil Literature: 1. Luurtsema G. et al. Fully automated high yield synthesis of (R)- and (S)-[¹¹ C]verapamil for measuring p-glycoprotein function with positron emission tomography. J. Labelled Compds. Radiopharm. 2000, 45, 1199-1207. 2. Wegmann T.D. et al. An improved method for the preparation of [¹¹ C]verapamil. Appl. Radiat. Isot. 2002, 57, 505-507.	5610.0000.5 0.5 mg per vial 5610.0001: 1 mg per vial 5610.0002: 2 mg per vial Please inquire for customized filling and bulk quantities. 

date of product catalogue issue: 05 April 2012