


Catalogue Number	Product	Order number / Unit
6160	Fluoroethylcholine chloride Reference standard for [¹⁸ F]Fluoroethylcholine Molar Mass: 171.64 C ₆ H ₁₅ ClFNO [479407-08-4] Yellowish solid packaged in dark glass screw cap vials. Purity: > 95 % Certificates: CoA; ¹ H and ¹⁹ F NMR spectra Chemical Name: CA index name: Ethanaminium, N-(2-fluoroethyl)-2-hydroxy-N,N-dimethyl-, chloride Synonymes: Fluoroethylcholine chloride; FECh chloride; FEtCh chloride Literature: 1. Hara T. et al. Development of ¹⁸ F-fluoroethylcholine for cancer imaging with PET: synthesis, biochemistry, and prostate cancer imaging J. Nucl. Med. 2002, 43, 187-199. 2. DeGrado T.R. et al. Synthesis and Evaluation of ¹⁸ F-Labeled Choline Analogs as Oncologic PET Tracers J. Nucl. Med. 2001, 42, 1805-1814.	6160.0010: 10 mg per vial Please inquire for customized filling and bulk quantities. 

date of product catalogue issue: 05 April 2012