

Catalogue  
Number

Product

Order number / Unit

---

**6176**

**Dibromomethane**

6176.0000.5: 0.5 ml per vial  
Please inquire for customized  
filling and bulk quantities.

**Molar Mass:** 173.83

CH<sub>2</sub>Br<sub>2</sub>

[74-95-3]

Colourless liquid packaged in dark glass crimp cap vials.

**Purity:** > 95 %

**Certificates:**

CoA; <sup>1</sup>H NMR spectrum

**Chemical Name:**

CA index name: Methane, dibromo-

**Synonymes:**

Dibromomethane; Methylene bromide; Methylene dibromide

**Literature:**

no literature reference available

