


Catalogue Number	Product	Order number / Unit
6179	1,2-Bis(nosyloxy)ethane Precursor for 1-[¹⁸F]Fluoro-2-nosyloxy-ethane Molar Mass: 432.38 $C_{14}H_{12}N_2O_{10}S_2$ [25297-82-9] Colourless to yellow crystals packaged in dark glass screw cap vials. Purity: > 95 % Certificates: CoA; ¹ H NMR spectrum Chemical Name: CA index name: Benzenesulfonic acid, 4-nitro-, 1,2-ethanediyl ester Synonymes: Benzenesulfonic acid, p-nitro-, ethylene ester; Ethylene glycol, bis(p-nitrobenzenesulfonate); Ethylene dinosylate; Ethylene glycol dinosylate; 1,2-Ethylene dinosylate; 1,2-Bis(nosyloxy)ethane; 1,2-Bis(4-nitrobenzenesulfonyloxy)ethane; 1,2-Bis(p-nitrobenzenesulfonyloxy)ethane; 1,2-Ethandiol dinosylate; Di-O-nosylglycol; Literature: 1. Musachio J.L. et al. Radiosyntheses and reactivities of novel [¹⁸ F]2-fluoroethyl arylsulfonates. J. Labelled Compd. Radiopharm. 2005, 48, 735-747.	6179.0100: 100 mg per vial 6179.1000: 1 g per vial Please inquire for customized filling and bulk quantities. 

date of product catalogue issue: 05 April 2012