


Catalogue Number	Product	Order number / Unit
6184	2-Fluoroethyl tosylate Reference standard for 2-[¹⁸F]Fluoroethyl tosylate Molar Mass: 218.25 $C_9H_{11}FO_3S$ [383-50-6] Colourless to pink liquid packaged in dark glass screw cap vials. Purity: > 95 % Certificates: CoA; ¹ H and ¹⁹ F NMR spectra Chemical Name: CA index name: Ethanol, 2-fluoro-, 1-(4-methylbenzenesulfonate) Synonymes: Ethanol, 2-fluoro-, 4-methylbenzenesulfonate; Ethanol, 2-fluoro-, p-toluenesulfonate; 1-Fluoro-2-(tosyloxy)ethane; 2-Fluoroethyl 4-methylbenzenesulfonate; 2-Fluoroethyl 4-toluenesulfonate; 2-Fluoroethyl p-toluenesulfonate; 2-Fluoroethyl tosylate; 4-Methylbenzenesulfonic acid 2-fluoroethyl ester; Fluoro(tosyloxy)ethane; p-Toluenesulfonic acid 2-fluoroethyl ester; Toluene-4-sulfonic acid 2-fluoroethyl ester Literature: 1. Musachio J.R. et al. Radiosyntheses and reactivities of novel [¹⁸ F]2-fluoroethyl arylsulfonates. J. Labelled Comp. Radiopharm. 2005, 48, 735-747.	6184.0100: 100 mg per vial 6184.1000: 1 g per vial Please inquire for customized filling and bulk quantities. 

date of product catalogue issue: 05 April 2012