

Catalogue Number	Product	Order number / Unit
6186	1,4-Bis(tosyloxy)but-2-yne Precursor for 1-[¹⁸ F]Fluoro-4-tosyloxy-but-2-yne Molar Mass: 394.46 C ₁₈ H ₁₈ O ₆ S ₂ [6337-59-3] Brown solid packaged in dark glass screw cap vials. Purity: > 95 % Certificates: CoA; ¹ H NMR spectrum Chemical Name: CA index name: 2-Butyne-1,4-diol, 1,4-bis(4-methylbenzenesulfonate); 2-Butyne-1,4-diol, bis(4-methylbenzenesulfonate); p-Toluenesulfonic acid, 2-butynylene ester Synonyms: 1,4-Bis(4-methylphenylsulphonyloxy)but-2-yne; 2-Butyne-1,4-diol ditosylate; 1,4-bis(tosyloxy)-but-2-yne; 1,4-Bis(tosyloxy)-2-butyne; 1,4-Butyne diol ditosylate; 1,4-Butynediol ditosylate; 1,4-Di-p-tosyl-2-butyne Literature: no literature reference available	6186.0100: 100 mg per vial 6186.1000: 1 g per vial Please inquire for customized filling and bulk quantities. 