


Catalogue Number	Product	Order number / Unit
6190	3-Fluoropropyltosylate Reference standard for 3-[¹⁸F] Fluoropropyltosylate Molar Mass: 232.27 $C_{10}H_{13}FO_3S$ [312-68-5] Colourless liquid packaged in dark glass screw cap vials. Purity: > 95 % Certificates: CoA; ¹ H NMR and ¹⁹ F NMR spectra Chemical Name: CA index name: 1-Propanol, 3-fluoro-, 1-(4-methylbenzenesulfonate); 1-Propanol, 3-fluoro-, 4-methylbenzenesulfonate; 1-Propanol, 3-fluoro-, p-toluenesulfonate Synonymes: 1-Fluoro-3-tosyloxypropane; 3-Fluoropropyl 4-toluenesulfonate; Toluene-4-sulfonic acid 3-fluoropropyl ester Literature: 1. Koivula T. et al. Simplified synthesis of N-(3-[¹⁸ F]fluoropropyl)-2β-carbomethoxy-3β-(4-fluorophenyl) nortropane ([¹⁸ F]β-CFT-FP) using [¹⁸ F]fluoropropyl tosylate as the labelling reagent. Labelled Compd. Radiopharm. 2005, 48, 463-471.	6190.0100: 100 mg per vial 6190.1000: 1 g per vial Please inquire for customized filling and bulk quantities. 

date of product catalogue issue: 05 April 2012