
Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

6200 Ethyl (p-tosyloxy)acetate
Precursor for [¹⁸F]Fluoroacetate

Molar Mass: 258.29

C₁₁H₁₄O₅S

[39794-75-7]

Colourless solid packaged in dark glass screw cap vials.

Purity: > 95 %

Certificates:

CoA; ¹H and ¹³C NMR spectra

Chemical Name:

CA index name: Acetic acid, [[(4-methylphenyl)-sulfonyl]oxy]-, ethyl ester


Synonymes:

Tosylethylglycolate

Literature:

1. Fujibayashi Y. et al. New approach to fully automated synthesis of sodium [¹⁸F]fluoroacetate - a simple and fast method using a commercial synthesizer. Nucl. Med. Biol. 2006, 33, 153-158.
2. Jeong J.M. et al. Synthesis of no-carrier-added [¹⁸F]fluoroacetate. J. Labelled Compd. Radiopharm. 1997, 39, 395-399.

6200.0010: 10 mg per vial
6200.0100: 100 mg per vial
Please inquire for customized filling and bulk quantities.



date of product catalogue issue: 05 April 2012