


Catalogue Number	Product	Order number / Unit
6210	Potassium fluoroacetate Reference standard for [¹⁸ F]Fluoroacetate Molar Mass: 116.13 C ₂ H ₂ FKO ₂ [23745-86-0] Colourless solid packaged in dark glass screw cap vials. Purity: > 95 % Certificates: CoA; ¹ H, ¹⁹ F, and ¹³ C NMR spectra Chemical Name: CA index name: Acetic acid, fluoro-, potassium salt Synonymes: Acetic acid, fluoro-, potassium salt Literature: 1. Fujibayashi Y. et al. New approach to fully automated synthesis of sodium [¹⁸ F]fluoroacetate - a simple and fast method using a commercial synthesizer. Nucl. Med. Biol. 2006, 33, 153-158. 2. Jeong J.M. et al. Synthesis of no-carrier-added [¹⁸ F]fluoroacetate. J. Labelled Compd. Radiopharm. 1997, 39, 395-399.	6210.0100: 100 mg per vial Please inquire for customized filling and bulk quantities. 

date of product catalogue issue: 05 April 2012