


Catalogue Number	Product	Order number / Unit
7100	S-Benzoyl-MAG-3 Ligand for labelling with ^{99m}Tc Molar Mass: 367.38 C₁₅H₁₇N₃O₆S [103725-47-9] Colourless crystals packaged in dark glass screw cap vials. Purity: > 95 % Certificates: CoA; ¹ H NMR spectrum Chemical Name: CA index name: Glycine, N-[(benzoylthio)acetyl]glycylglycyl- Synonymes: Glycine, N-[N-[N-[(benzoylthio)acetyl]glycyl]glycyl]; Betiatide; N-(S-Benzoylmercapto)acetylglycylglycylglycine; S-Benzoylmercaptoacetyltriglycine; MP 600; S-Benzoyl-MAG-3; BzSAC-Gly-Gly-Gly-OH	7100.0100: 100 mg per vial 7100.1000: 1 g per vial Please inquire for customized filling and bulk quantities. 
	Literature: 1. Bormans G. et al. Investigation of the labeling characteristics of ^{99m} Tc- mercaptoacetyltriglycine. Nucl. Med. Biol. 1995, 22, 339-349. 2. Noll B. et al. Preparation of the renal function and imaging agent technetium-99m-MAG3 starting from S-unprotected mercaptoacetyltriglycine. Appl. Radiat. Isot. 1992, 43, 899-901.	

date of product catalogue issue: 05 April 2012