
Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

7150**BMEDA****Ligand for labelling with ^{99m}Tc****Molar Mass:** 236.44C₁₀H₂₄N₂S₂

[93798-62-0]

Colourless to yellowish liquid packaged in dark glass crimp cap vials.

Purity: > 90 %**Certificates:**CoA; ¹H NMR spectrum; GC**Chemical Name:**


CA index name: Ethanethiol, 2,2'-[[2-(diethylamino)ethyl]imino]bis-

Synonymes:

N,N-Bis(2-mercaptoethyl)-N',N'-diethylethylenediamine

Literature:

1. Bao A. et al. Direct ^{99m}Tc labelling of pegylated liposomal doxorubicin (doxil) for pharmacokinetic and non-invasive imaging studies. J. Pharmacol. Exp. Ther. 2004, 308, 419 -25.
2. Leon A. et al. Novel mixed ligand technetium complexes as 5-HT1A receptor imaging agents. Nucl. Med. Biol. 2002, 29, 217-26.

7150.0010: 10 mg per vial
7150.0100: 100 mg per vial
Please inquire for customized filling and bulk quantities.

date of product catalogue issue: 05 April 2012