

Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

**7165 EDTMP****Ligand for radiolabelling****Molar Mass:** 436.12 $C_6H_{20}N_2O_{12}P_4 \cdot x H_2O$ 

[1429-50-1]

White to off-white powder packaged in dark glass screw cap vials.

**Purity:** > 95 %**Certificates:**CoA; <sup>1</sup>H spectrum**Chemical Name:**

CA index name: Phosphonic acid, P,P',P'',P'''-[1,2-ethanediylbis[nitrilobis(methylene)]]tetrakis-

**Synonyms:**

Phosphonic acid, [1,2-ethanediylbis[nitrilobis(methylene)]]tetrakis-; Phosphonic acid, [ethylenebis(nitrilodimethylene)]tetra-; Briquest 422-25S; Chelating agent 540; Chelest PH 540; Cublen 115; Dequest 2040; Dequest 2041; Editempa; EDPA; EDPA (chelating agent); EDTF; EDTMP; EDTMPA; EDTPA; EDTPH; Ethylenedi(nitrilodimethylene)tetrakisphosphonic acid; Ethylenediamine-N,N,N',N'-tetra(methylphosphonic acid); Ethylenediamine-N,N,N',N'-tetrakis(methylenephosphonic acid); Ethylenediaminetetra(methylenephosphonic acid); Ethylenediaminetetrakis(methylenephosphonic acid); Ethylenediaminetetrakis(methylphosphonic acid); Ethylenediaminetetramethylenephosphonate; Ethylenediaminetetra(methylenephosphonic acid); N,N,N',N'-Tetrakis(phosphonomethyl)ethylenediamine

**Literature:**

no literature reference available

Please inquire for customized filling and bulk quantities.

