

Catalogue Number	Product	Order number / Unit
7211	<p><b>Copper tetraMIBI tetrafluoroborate (GMP)</b></p> <p><b>Ligand for labelling with <math>^{99m}\text{Tc}</math> or <math>^{188}\text{Re}</math></b></p> <p>Manufactured according to GMP requirements for APIs (ICH Q7) European Drug Master File (EDMF) available.</p> <p><b>Molar Mass:</b> 602.98 <math>\text{C}_{24}\text{H}_{44}\text{CuN}_4\text{O}_4 \cdot \text{BF}_4</math> [103694-84-4]</p> <p>Colourless to white crystalline solid packaged in dark glass screw cap vials.</p> <p><b>Purity:</b> &gt; 99 %</p> <p><b>Certificates:</b> CoA with <math>^1\text{H}</math>, <math>^{13}\text{C}</math>, and <math>^{19}\text{F}</math> NMR and IR spectra and melting point (identity); HPLC (purity); HS-GC (residual solvents); Assay (HPLC); Assay Boron and Copper (ICP-MS); tested for Bioburden and Bacterial Endotoxins</p> <p><b>Chemical Name:</b> For chemical information please refer to product number 7210.</p> <p><b>Synonymes:</b> Please refer to product number 7210.</p> <p><b>Literature:</b> For literature references please refer to product number 7210.</p>	<p>7211.0100: 100 mg per vial 7211.1000: 1 g per vial Please inquire for customized filling and bulk quantities.</p> 