

Catalogue Number	Product	Order number / Unit
802	<p>Lithium aluminium hydride 0.1 M 0.1 M solution in dried tetrahydrofuran with controlled ¹²C-background</p> <p>Manufactured and packaged in argon atmosphere. Repeated extractions from vial in argon atmosphere does not affect product quality.</p> <p>Molar Mass: 37.95</p> <p>LiAlH₄</p> <p>[16853-85-3]</p> <p>Clear colourless liquid packaged in clear glass vials (2 ml) sealed with teflon-faced rubber stoppers and tear-off crimp caps.</p> <p>Purity: < 5 ppm background CO₂ 50 ppm Methanol</p> <p>Certificates: CoA; Background ¹²C ([¹²C]Carbonate determined as [¹²C]CO₂ and [¹²C]Methanol as methyl acetate by gas chromatography)</p> <p>Chemical Name: Please refer to product number 801.</p> <p>Synonymes: Please refer to product number 801.</p> <p>Literature: Please refer to product number 801.</p>	<p>802.0001: 1 ml per vial Please inquire for customized filling and bulk quantities.</p> <p>LiAlH₄</p>