

Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

**808 TBA HCO<sub>3</sub> 0.075 M**

**Aqueous solution, stabilized with ethanol  
Tested for sterility and bacterial endotoxins  
according to Ph. Eur.**

Manufactured according to GMP requirements for APIs (ICH Q7)

**Molar Mass:** 303.48

C<sub>17</sub>H<sub>37</sub>NO<sub>3</sub>

[17351-62-1]

Clear colourless liquid packaged in PP vials with crimp cap.

**Purity:**

**Certificates:**

CoA; pH; bacterial endotoxines/LAL test; test for sterility; Assay

**Chemical Name:**

CA index name: Tetrabutylammonium bicarbonate

**Synonymes:**

Tetrabutylammonium hydrogen carbonate solution

**Literature:**

TBA solution is designed as an reagent for the synthesis of [<sup>18</sup>F]FDG using the GE TRACERlab FXFDG Synthesizer(r) module (alkaline or acidic TBA method).

**808.0000.6: 0.6 ml per vial  
Please inquire for customized  
filling and bulk quantities.**

