

Catalogue Number	Product	Order number / Unit
830	Hydriodic acid Redistilled Aqueous solution, > 55 %, redistilled Manufactured and packaged in argon atmosphere. Storage at 5±3 degC Molar Mass: 127.91 HI [10034-85-2] Clear liquid, colourless to slightly yellow packaged in amber headspace glass vials (10 ml) sealed with teflon-coated rubber stoppers and tear-off crimp caps. Purity: > 99.9 wt % Certificates: CoA; Assay (titration) Chemical Name: CA index name: Hydriodic acid Synonymes: Hydriodic acid; Hydrogen iodide; Hydrogen monoiodide Literature: no literature reference available	830.0003.5: 3.5 ml per vial 830.0007: 7 ml per vial Please inquire for customized filling and bulk quantities. HI