

Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

9602 Scrambled Ubiquicidin (29-41)

9602.0001: 1 mg per vial
9602.0010: 10 mg per vial
Please inquire for customized filling and bulk quantities.

Molar Mass: 1692.95

$C_{68}H_{121}N_{31}O_{18}S$

CAS-RN not yet assigned

Colourless to off-white freeze-dried solid packaged in dark glass screw cap vials.

Purity: $\geq 95\%$

Certificates:

CoA; MS (identity) and HPLC (purity)

Chemical Name:


L-Threonine, L-lysyl-L-arginyl-L-asparaginyl-L-glutaminyl-L-arginyl-L-methionyl-L-alanyl-L-arginyl-L-tyrosyl-L-arginyl-L-arginylglycyl-

Synonymes:

(29-41)-Ubiquicidin, scrambled

Literature:

1. Welling M.M. et al. Radiochemical and biological characteristics of ^{99m}Tc -UBI 29-41 for imaging of bacterial infections. Nucl. Med. Biol. 2002, 29, 413-422.



9602 scrambled Ubiquicidin (29-41)
H-Lys-Arg-Ala-Gln-Arg-Met-Ala-Arg-Tyr-Arg-Arg-Gly-Thr-OH

date of product catalogue issue: 05 April 2012