

Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

**9730 HYNIC-NOC trifluoroacetate**

Precursor for [<sup>99m</sup>Tc]HYNIC-NOC

**Molar Mass:** 1204.42

C<sub>59</sub>H<sub>73</sub>N<sub>13</sub>O<sub>11</sub>S<sub>2</sub>

CAS-RN not yet assigned

Colourless freeze-dried solid packaged in dark glass screw cap vials.

**Purity:** ≥ 95 %

**Certificates:**

CoA; MS (identity) and HPLC (purity)

**Chemical Name:**

L-Cysteinamide, N-[(6-hydrazinyl-3-pyridinyl)carbonyl]-D-phenylalanyl-L-cysteinyl-L-(1-naphthyl)alanyl-D-tryptophyl-L-lysyl-L-threonyl- N-[(1R,2R)-2-hydroxy-1-(hydroxymethyl)propyl]-, cyclic (2→7)-disulfide


**Synonymes:**

L-Cysteinamide, N-[(6-hydrazino-3-pyridinyl)carbonyl]-D-phenylalanyl-L-cysteinyl-L-(1-naphthyl)alanyl-D-tryptophyl-L-lysyl-L-threonyl-N-[(1R,2R)-2-hydroxy-1-(hydroxymethyl)propyl]-, cyclic (2→7)-disulfide; HYNIC-[1-Nal<sup>3</sup>] Octreotide

**Literature:**

no literature reference available

9730.0001: 1 mg per vial  
 9730.0010: 10 mg per vial  
 Please inquire for customized filling and bulk quantities.



HYNIC-NOC trifluoroacetate  
 HYNIC-o-Phe-Cys-1-Nal-o-Trp-Lys-Thr-Cys-Thr(o)