

Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

9762 **NOTA-Octreotide trifluoroacetate**
Precursor for radiolabeled NOTA-Octreotide

Molar Mass: 1304.50

$C_{61}H_{85}N_{13}O_{15}S_2 \cdot x C_2HF_3O_2$

CAS-RN not yet assigned

Colourless to off-white freeze-dried solid packaged in dark glass screw cap vials.

Purity: $\geq 95\%$

Certificates:

CoA; MS (identity) and HPLC (purity)

Chemical Name:

NOTA-D-Phe-cyclo[Cys-Phe-DTrp-Lys-Thr-Cys]-Throl


Synonymes:

NOTA-octreotide; IMP466

Literature:

1. Laverman P. et al. A Novel Facile Method of Labeling Octreotide with ^{18}F -Fluorine. J. Nucl. Med. 2010, 51, 454-461.

9762.0001: 1 mg per vial
9762.0010: 10 mg per vial
Please inquire for customized filling and bulk quantities.



NOTA-Octreotide trifluoroacetate

NOTA-o-Phe-Cys-Phe-o-Trp-Lys-Thr-Cys-Thr(ol)