

Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

**9804 NODAGA-RGD dimer acetate****Precursor for radiolabeled NODAGA-RGD dimer****Molar Mass:** 1707.80 $C_{74}H_{110}N_{22}O_{25}$ 

CAS-RN not yet assigned

Colourless to off-white freeze-dried solid packaged in dark glass screw cap vials.

**Purity:**  $\geq 95\%$ **Certificates:**


CoA; MS (identity) and HPLC (purity)

**Chemical Name:**Cyclo(L-arginylglycyl-L- $\alpha$ -aspartyl-D-tyrosyl-L-lysyl), N6-([4,7-bis(carboxymethyl)octahydro-1H-1,4,7-triazonin-1-yl]acetyl)-L-glutamoyl]bis-**Synonymes:**

NODAGA-RGD

**Literature:**

no literature reference available

**9804.0001: 1 mg per vial**  
**9804.0010: 10 mg per vial**  
Please inquire for customized filling and bulk quantities.

date of product catalogue issue: 05 April 2012