


Catalogue Number	Product	Order number / Unit
9805	<p>NODAGA-RGD trifluoroacetate</p> <p>Precursor for radiolabelled RGD peptides</p> <p>Molar Mass: 977.03</p> <p>$C_{42}H_{64}N_{12}O_{15}$</p> <p>CAS-RN not yet assigned</p> <p>Colourless to off-white solid packaged in dark glass screw cap vials.</p> <p>Purity: $\geq 95\%$</p> <p>Certificates: CoA; MS (identity) and HPLC (purity)</p> <p>Chemical Name: Cyclo[L-arginylglycyl-L-α-aspartyl-D-tyrosyl-N6-([[4,7-bis(carboxymethyl)octahydro-1H-1,4,7-triazonin-1-yl]-4-carboxy-butryl])-L-lysyl]</p> <p>Synonymes: Cyclo[L-arginylglycyl-L-α-aspartyl-D-tyrosyl-N6-([[4,7-bis(carboxymethyl)-1,4,7-triazonan-1-yl]acetyl])-L-lysyl]; NODAGA-RGD</p> <p>Literature: no literature reference available</p>	<p>9805.0001: 1 mg per vial 9805.0010: 10 mg per vial Please inquire for customized filling and bulk quantities.</p>  <p>NODAGA-RGD trifluoroacetate</p> <p>Arg-Gly-Asp-D-Tyr-Lys(NODAGA)</p>