

Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

9836 HYNIC-Lanreotide trifluoroacetate
Precursor for radiolabelled HYNIC-Lanreotide

Molar Mass: 1231.45

$C_{60}H_{74}N_{14}O_{11}S_2$

CAS-RN not yet assigned

Colourless to off-white freeze-dried solid packaged in dark glass screw cap vials.

Purity: $\geq 95\%$

Certificates:

CoA; MS (identity) and HPLC (purity)

Chemical Name:

L-Threoninamide, 3-(2-naphthalenyl)-N-[(6-hydrazinyl-3-pyridinyl)carbonyl]-D-alanyl-L-cysteinyl-L-Tyrosyl-D-tryptophyl-L-lysyl-L-valyl-L-cysteinyl-, cyclic (2→7) disulfide


Synonymes:

HYNIC-Lanreotide

Literature:

no literature reference available

9836.0001: 1 mg per vial
9836.0010: 10 mg per vial
Please inquire for customized filling and bulk quantities.



HYNIC-Lanreotide trifluoroacetate
HYNIC-D-Phe-Cys-Tyr-D-Trp-Lys(Boc)-Thr-Cys-Thr-NH₂