

Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

9860 **Cyclo-RGDfK**
Integrin targeting sequence.

Molar Mass: 603.67

$C_{27}H_{42}N_9O_7$

[161552-03-0]

Colourless freeze-dried solid packaged in dark glass screw cap vials.

Purity: $\geq 95\%$

Certificates:

CoA; MS (identity) and HPLC (purity)

Chemical Name:

CA index name: Cyclo(L-arginylglycyl-L- α -aspartyl-D-phenylalanyl-L-lysyl)


Synonymes:

cyclo(-Arg-Gly-Asp-D-Phe-Lys-)

Literature:

1. Liu S. Radiolabeled Cyclic RGD Peptides as Integrin $\alpha v\beta 3$ Targeted Radiotracers: Maximizing Binding Affinity via Bivalency. *Bioconjugate Chem.*, 2009, 20, 2199-2213.

9860.0001: 1 mg per vial
9860.0010: 10 mg per vial
Please inquire for customized filling and bulk quantities.



Cyclo-RGDfK

— Arg-Gly-Asp-D-Phe-Lys —