

Catalogue Number	Product	Order number / Unit
------------------	---------	---------------------

9903 [NaI³,Lys⁵(Boc)]octreotide acetate

Precursor for N-terminal substituted [NaI³] octreotide analogs

Molar Mass: 1169.42

C₅₈H₇₆N₁₀O₁₂S₂

CAS-RN not yet assigned

Colourless to off-white freeze-dried solid packaged in dark glass screw cap vials.

Purity: ≥ 95 %

Certificates:

CoA; MS (identity) and HPLC (purity)

Chemical Name:

L-Cysteinamide, D-phenylalanyl-L-cysteinyl-3-(1-naphthalenyl)-L-alanyl-D-tryptophyl-N6-[(1,1-dimethylethoxy)carbonyl]-L-lysyl-L-threonyl-N-[(1R,2R)-2-hydroxy-1-(hydroxymethyl)propyl]-, cyclic (2→7)-disulfide


Synonymes:

[NaI³,Lys⁵(Boc)]octreotide

Literature:

no literature reference available

9903.0001: 1 mg per vial
 9903.0010: 10 mg per vial
 Please inquire for customized filling and bulk quantities.



[Lys⁵(Boc)]lanreotide acetate

H-D-2-Nal-Cys-Tyr-D-Trp-Lys(Boc)-Val-Cys-Thr-NH₂